

TECHNICAL DATA

# External vibration generator for PHYWE Ripple Tank

Article no: 11260-10



## Function and applications

Optional accessory for ripple tank: second dipper for demonstrating interference patterns from waves that are not in phase.

## Benefits

- Settings and power for this dipper come directly from the ripple tank
- No additional power supply needed
- Stand for the second dipper (padded foam base) is included in the set

## Equipment and technical data

- External vibration generator (dipper) for ripple tank on damped base

## Scope of delivery

Connecting cable, 32 A, black, various lengths	07361-05	2
--	----------	---